

A Novel Epigenomic Profiling Assay Reveals Disease Biology and Therapeutic Targets in Gastro-esophageal Cancer from 1mL of Plasma

Jonathan Beagan,¹ Aparna Gorthi,¹ Anthony D'Ippolito,¹ Travis Clark,¹ Michael Coyne,¹ Tyrone Tamakloe,¹ Amy Donahue,¹ Baovy Tran,¹ Hyun-Hwan Jeong,¹ Kristian Cibulskis,¹ Hathairat Sawaengsri,¹ Illexa Ashley Schechter², Corrie A. Painter,¹ Juliann Chmielecki,¹ J Carl Barrett¹ Mosammat Faria Afreen,² Aditya Pandey,² Martin S. Taylor,² Matthew R. Strickland,² Matthew Eaton,² Samuel J. Klempner²

¹Precede Biosciences, Boston, MA ²Mass General Cancer Center, Boston, MA

ABSTRACT # 8833

INTRODUCTION

- Gastro-esophageal adenocarcinoma (GEA) is an aggressive and heterogeneous disease with a poor prognosis.
- Targeted therapies have failed largely due to challenges associated with biopsy collection, tumor antigen heterogeneity, and accurately distinguishing transcriptionally distinct esophageal, GE junction or gastric tumor types.
- We developed a novel, multimodal liquid biopsy assay that profiles genome-wide transcriptional activation¹.
- Here, we deploy this assay to define clinically relevant insights in a cohort of GEA patients.

METHODS


We profiled 107 1mL plasma samples from 65 patients with gastroesophageal cancer enrolled in an IRB-approved protocol Massachusetts General Hospital (MGH). 56 samples represented pre-treatment blood draws with matched HER2 IHC pathologist assessment. 46 pre-treatment samples yielded high-quality datasets across all analytes (promoter, enhancer, DNAm). 15 samples had detectable ctDNA by ichorCNA estimation². HER2 IHC, staging, and histology grade breakdown of this cohort is below.

HER2 IHC (all samples, n=46/ detectable ctDNA by ichor, n=15)			
0	1+	2+	3+
26 / 7	5 / 2	8 / 2	7 / 4
Staging at Sample Collection (all samples, n=46/ detectable ctDNA by ichor, n= 15)			
III		IV	
27 / 3		19 / 12	
Histology grade (all samples, n=46/ detectable ctDNA by ichor, n=15)			
Poorly Differentiated or Grade 3	Moderate-poorly Differentiated or Grade 2	Not Reported	
12 / 4	11 / 2	23 / 9	

We developed a HER2 classifier trained on epigenomic features (enhancers, promoters, DNAm peaks) which enabled classification of HER2 ± cancer cell lines, as previously presented (Parsons et al., SABCS 2023)³.


- We further tuned this model for application to gastroesophageal cancer plasma samples.
- Classifier was cross-validated using standard leave-one-out scheme.

Comprehensive epigenomic platform offers dynamic resolution into target and pathway biology from 1 mL of plasma



Cell free DNA derived from tumors exist in circulation as chromatin fragments that faithfully maintain tumor-associated epigenetic modifications on the histones and DNA. Antibodies against H3K27ac marking active enhancers, H3K4me3 marking active promoters and DNA methylation are used to enrich for associated DNA fragments from 1 mL plasma and sequenced to define genome-wide epigenomic maps that capture the underlying transcriptional state of the tumor cells.¹


Figure 1: Plasma-based epigenomic activation signal at GEA marker genes



Genome-wide promoter and enhancer regions were assessed for correlation of signal with ichor ctDNA% estimates. **Top left:** Volcano plot showing genes with positive fold change display increased signal in GEA cancer samples (red). **Top Right, Bottom:** Promoter tracks and scatterplots of promoter/enhancer signal vs estimated ctDNA% (above) are displayed for select growth factors (*FGFR2*, *VEGFA*) and transcription factors (*CDX2*, *EVX1*, *SOX9*) previously associated with gastroesophageal cancer.

RESULTS

Figure 2: Epigenomic activation panel to guide patient selection for ADC targets



Left: Patient-specific (by row) promoter signals at genes representing select ADC targets (plus housekeeping genes *GAPDH/ACTB*). Samples are sorted by ichor ctDNA fraction estimate. **Middle:** Enhancer and promoter signals for each gene were ctDNA-corrected and combined into an epigenomic activation score, represented on a white (low) to red (high) gradient. **Right:** LBx detection of elevated *MET* signal initially absent from tissue biopsy of Patient A. Elevated *MET* promoter signal identified by promoter tracks and promoter signal vs ctDNA linear residual (scatter plot). Patient B is a high ctDNA GEA patient sample with moderate-low *MET* promoter signal included for comparison.

Figure 3: Evaluating epigenomic status of modulators of ADC response







Figure 4: Multi-analyte HER2 classifier robustly stratifies patients



Top Left: ROC curve assessing performance of regularized logistic regression model to classify HER2 IHC status (positive vs negative)³ of samples with detectable ctDNA by ichor. **Top Right:** Stratification of the model's output probabilities by HER2 IHC reveal one IHC 0 patient assigned a high likelihood of being HER2 positive. **Bottom:** Performance of HER2 classifier, assessed via AUC, for all samples and those with detectable ctDNA after removal of outlier IHC 0 patient for which a HER2 positive classification may have been beneficial.

CONCLUSIONS

- Here, we present a novel liquid biopsy approach that leverages comprehensive epigenomic profiling from 1 mL of plasma, an approach that represents an advance over current tissue and liquid based methods.
- Our genome-wide approach enabled direct comparison of the activation status of GEA-relevant ADC target genes such as *HER2*, *HER3*, and *MET* along with genes indicating sensitivity to warhead class, highlighting potential for use in patient selection.
- Using a HER2 classifier trained on enhancer, promoter, and DNA methylation profiles, we accurately stratify patients based on IHC/FISH status, identifying those who may benefit from *HER2*-directed therapy.
- Our minimally invasive approach enables longitudinal expression profiling of therapy related dynamics, thereby informing rational sequential therapeutic strategies.